

hangzhou onicon kio-lab center with most advanced kio-laboratory, pilot plant and producing facilities as well as all necessary analytical equipments(NMR,HPLC, GS,AXIMA CFR PLUS,ETC) is willing to offer the below working standard sample and test samples of API.

Alarelin Acetate	cas	79561-22-1	Atosiban Acetate	cas	90779-69-4
GHRP-6 Acetate	cas	87616-84-0	Melanotan II Acetate	cas	121062-08-6
Octreotide Acetate	cas	83150-76-9	Calcitonin(Salmon)	cas	47931-85-1
GHRP-2 Acetate	cas	158861-67-7	Bivalirudin Acetate	cas	128270-60-0
Leuprorelin Acetate	cas	53714-56-0	Eptifibatide Acetate	cas	148031-34-9
Cetrorenlin Acetate	cas	130143-01-0	Tetracosactide Acetate(ACTH1-24)	cas	16960-16-0
Triptorelin Acetate	cas	57773-63-4	Desmopressin Acetate	cas	16789-98-3
Exenatide Acetate	cas	141732-76-5	Gonadorelin Acetate	cas	34973-08-5
Nesiritide Acetate	cas	114471-18-0	Sermorelin Acetate	cas	86168-78-7
Pramlintide Acetate	cas	196078-30-5	Terlipressin Acetate	cas	14636-12-5
Argipressin Acetate	cas	113-79-1	Thymopentin	cas	69558-55-0
Lysipressin Acetate	cas	50-57-7	PT-141 Acetate	cas	32780-32-8
Teriparatide Acetate	cas	52232-67-4	Deslorelin Acetate	cas	57773-65-6

Hangzhou Onicon Chemical co.,ltd	
ADD:	3F Hangzhou Pharmaceutical(Yiyao) Building, No.439 Zhong Shan North Rd., Hangzhou China
TEL:	+86-571-85054016
FAX:	+86-571-85054082
URL:	Http://www.onicon.cn
EMAIL:	cjw9999@163.com